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Waste Analysis using Ion
Chromatography Method

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Presently, criteria and practise to control potential harms of waste materials to the environment are defined by the Ministry of Environment and Urbanization in a Legislation Concerning Waste. Common purpose of waste regulations can be defined as prohibition of all activities including direct or indirect disposal, storage, transfer, removal etc. of all types of waste and effluent in a manner to give potential hazard or harm to the environment, control of consumer goods which could be harmful for the environment and to save the animal and plant generations as well as to keep ecological balance by preventing hazardous contaminants and pollutants for air, water and soil. Waste regulations include defining, applying and inventing principles, policies and programmes in this scope.

When Annex 1 of Tests and Sampling Methods for solid waste analysis and Annex 2 for the Waste Acceptance Criteria as well as Annex 1 for Waste Oil Categories and Limit Values for Allowed Contaminant Parameters are examined, it is noticed that some toxic elements, ions and dangerous organic material have been restricted. Nowadays, ion chromatography technique is used to analyse Cl-, F- and SO42- anions for which limit values are determined by the respective legislation. These anions are known to cause toxic effects when mixed directly or indirectly after chemical modification with soil or water within ecological system.

Ion Chromatography is a chromatographic technique with a wide application area which allows the analysis of ions at ppm and ppb levels. This technique allows simultaneous, rapid and reliable analysis of F-, Cl-, SO42-, NO2-, NO3-, Br-, PO43-, ClO2-, ClO3-, BrO3- anions and Li+, Na+, Ca2+, K+, Mg2+, NH4+, cations in clean and waste water samples and all types of solid and liquid waste. Also known as ion exchange chromatography, for this technique, generally anion or cations exchange ion chromatography columns with packing materials of sulphonic acid or quarternary ammonium groups containing bound resin are used. Anions or cations to be analysed are separated in these columns with respect to their ionic charge, diameter of hydrate ions, polarization status, pH and concentration of the mobile phase and column temperature and analysed using conductivity detector. Although, in general highly sensitive conductivity detectors are used in ion chromatography instruments; for some specific analysis including cyanide, cyanogen

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
chloride and organic acids post-column derivatization or UV Absorbance detectors could also be used.

Below, is data from an application note prepared by Ant Teknik application specialists who used Shimadzu Prominence NS model ion chromatography system to analyse anions at limit values mentioned in Annex 2 of Waste Acceptance Criteria Legislation. Plastic solid waste samples are analysed using ion chromatography after the sample is burnt using calorimeter bomb. Table 1 shows chromatographic conditions for F-, Cl-, SO42- ions, Figure 1 shows graphic of detector signal versus concentration using calibration curve method and Figure 2 shows sample chromatogram for F-, Cl-, SO42- ions for 1 g of plastic solid waste sample.


Table 1. Chromatographic conditions for F-, Cl-, SO₄²⁻ ions

Instrument	: HPLC
Columns	: Shim-pack IC-A3 (150 mm x 4.6mm)
Column Oven Temp	: 40 °C
Injection Volume	: 50 uL
Flow Rate	: 1,2ml/min
Detector	: CDD
Time	: 20 minutes
Mobile Phase	: 3.2mM Bis-Tris 50mM Boric Acid

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Figure 1. External standard – graphics for detector signal versus concentration using calibration curve method

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Figure 2. Sample chromatogram for F-, Cl-, SO42- ions for 1 g of plastic solid waste sample

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References

1. Turkish Ministry of Environment and Urbanization, Legislation Concerning Solid Waste
2. EPA Method 9056A and EPA Method 5050
3. Ant Teknik - Shimadzu Application Note L006

About Ant Teknik

Ant Teknik, founded in 1999 in Istanbul, is offering analytical instruments, technical service and spare parts; validation, application and training as well turn-key laboratory project services for Quality Control and Research&Development laboratories. Today, we offer innovative analytical solutions for various fields such as environment, pharmaceutical, food and chemicals. We have offices in Istanbul, Ankara, Izmir, Adana and Baku; more than 90 employees; and regional subdealers.

We represent respected manufacturers of analytical and other laboratory instruments including Shimadzu Corporation, Kratos, Rudolph Research, JeioTech, Hitachi, Restek, GL Sciences in Turkey.

For detailed information about Waste Management and Environment Analysis Solutions from Ant Teknik, please visit www.antteknik.com.