Zafer Bayram

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*Recently, it was determined that colorants used in leather, textiles and garment manufacturing industry had adverse effects on human health. Azo dyes therefore, were taken to restricted chemicals list in European countries. In Turkey, starting from March 1st 1995, the usage of aryl amines and azo dyes containing aryl amines in leather and textiles processing was prohibited by the Ministry of Health. On the other hand, it is found out that some textile products especially imported from Far East contain banned colorants and controls and inspections by the Ministry are continuously done in order to prevent these.*


*Analytical instruments are used in textile laboratories for several purposes like the determination of banned colorants and the sources of problems on fabric materials as well as competitor product analyses or alternative product trials. Some basic applications done by HPLC, FTIR and AAS are shown below:*

### *1.p-Chloraniline Analysis by HPLC*

*Banned material analysis is one of the most important analysis performed on azo colorants. Dyed textile products should not contain p-chloraniline (PCA) over 20 ppm due to their carcinogenic effects. P-chloraniline analysis on dyed fabric sample is done by HPLC.*

*HPLC analysis of 30 ppm of PCA solution is shown on Figure 1.1. In this study, the characteristic structure of the PCA peak and wavelengths are shown. PCA peak is determined at around the 9th minute (retention time).*

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
Figure 1.1. HPLC analysis of 30 ppm of p-Chloroaniline standard

In order to perform PCA analysis by HPLC on dyed fabric sample, the colorant on the fabric material should be dissolved and transferred on a solvent media. PCA analysis is performed on this solvent media. Figure 1.2. demonstrates HPLC analysis on dyed fabric sample. This analysis is performed on the solvent obtained from a sample fabric dyed with an azo colorant, and a peak is detected with the same characteristic structure, wavelength and retention time. Therefore, it is concluded that the fabric sample dyed with this colorant contains PCA. PCA concentration is calculated using formula mentioned in test methods.

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Figure 1.2. HPLC analysis of colorant sample on dyed fabric material

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
### 2. Stain Analysis on Fabric by FTIR


FTIR Spectrophotometer is used to perform the analysis of stains on fabric cloth. The source of the stain occurred during manufacturing by FTIR analysis. Figure 2.1. shows FTIR spectra from different parts of a fabric cloth with and without stains. Looking at the interferograms from 2800-3000  $\text{cm}^{-1}$  wavelength range, a sharp peak is observed. This stain originated peak matches with the FTIR spectrum of lubricating oil shown on Figure 2.2. Therefore, it is concluded that this stain comes from an outer source (grease) from one of the machinery where the fabric cloth was previously processed on.

Other application areas for FTIR are competitor product analyses or alternative product trials.

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Figure 2.1. Comparison between FTIR spectra from different parts of fabric cloth with and without stains.

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Figure 2.2. FTIR spectrum of lubricating oil

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
### 3. Metal Analysis by Atomic Absorption Spectrophotometer (AAS)

Atomic absorption spectrophotometer is used to perform metal analyses on dyes, dyed fabric and unprocessed fabric samples. For the metal analysis on dye samples, the dye is diluted using the right solvent and the analysis is performed. For dyed and unprocessed fabric samples on the other hand, the samples should first be extracted in order to transfer the metal components to solvent media. Figure 3.1 shows the results of chromium analysis on fabric sample dyed with Reagent Black 8 dye sample. The analysis is done in comparison between fabric sample dyed with standard dye versus fabric sample dyed with sample dye.

Absorbance vs. concentration graphics should be drawn first for metal analyses. For this purpose, chromium samples with known concentration should be analysed by AAS and corresponding absorbance values should be determined. The chromium concentration of the samples is calculated using the obtained absorbance values and absorbance vs. concentration graphic equation. Real values are obtained by multiplying the result by the dilution factor.

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
Figure 3.1. Chromium analysis on fabric samples

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In AAS analyses, some differences might occur due to some factors like metal type and dye type. In iron analysis for example, standard samples were prepared as 0.5 ppm, 1 ppm and 2 ppm as shown in Figure 3.2.

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Figure 3.2. Iron analysis on fabric sample.

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