Serhat Çeşme

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Sales Engineer

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In the cosmetic industry, different analyzes are performed for various purposes such as quality control, product development, R&D studies and compliance with legal requirements. Perfume, cream and other cosmetic products as well as raw materials are analyzed using various analytical techniques and instruments.

Some examples are:

- Flavour and fragrance analysis: GCMS
- Preservatives (Paraben, Benzoic and Sorbic Acid, Triclosan, Formaldehyde etc): HPLC, GC, GCMS
- Hyaluronic acid: LCMSMS, SEC
- SPF (Sun Protection Factor): UV-VIS-NIR
- Lead: AAS
- Viscosity: Viscosimeter
- Density: Densitometer
- Concentration and Solids: Refractometer

In this article, viscosity, density and concentration measurements which are among the quality control parameters commonly used in the cosmetic industry, instruments used for these measurements and their other application areas will be explained.

Viscosity Measurements in Cosmetics

Viscosity is defined as the resistance of viscous materials against flow.

Viscosity measurements in the cosmetic industry are an important part of quality control both during the production process and for the finished products.

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For example, lip balm needs to be highly viscous for proper adhesion and protection of the skin. On the other hand, a body spray must have a very low viscosity enabling the product flow freely from the spray dispenser and show enhanced absorption through the skin.

The viscosity of a fluid is a parameter that must be measured while designing the pipeline or during the selection of the process equipments such as pumps and production line parameters such as pipe volume/diameter etc. For this reason, viscosity measurements have become a key factor for high quality production of cosmetics.

Fungilab V-Compact Viskozimetre

Fungilab, one of the world's leading viscometer manufacturers, is producing the award-winning V-Compact series, which is noted as "the world's first wireless access rotary viscometer". The instrument allows users to access test data and multiple application audits in different searches. The V-Compact enhances users' productivity and efficiency with features such as remote access via mobile application, Wi-Fi and Bluetooth functionality, and unique digital interface. It also provides remote online support in cases where urgent help is needed.

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RUDOLPH RESEARCH
DDM-2910 Model Densitometre

Application Areas of Viscosimeters:

- Cosmetic Industry (Cream, lotion, shampoo, detergent, nail polish etc.),
- Chemical Industry (Polymers, poliuretane based samples etc.),
- Paint Industry (paint, ceramic dyes) and Printing Industry (inks),
- Food Industry (Milk, yoghurt, chocolate, tomato paste, yeast etc.),
- Petrochemicals Industry (Mineral oil and derivatives, asphalt etc.),
- Pharmaceuticals Industry (Quality control testing of suspensions and emulsions, finishing etc.).

Density Measurements in Cosmetics and Fragrances

In flavour and fragrances industry, density is used to check the blended and/or diluted oils to the desired concentration and customer specifications. Raw materials are checked as well as alcohol, sugar (°Brix), % solids, and many other parameters. Specific Gravity is always checked on final product as a means of ensuring product consistency.

Rudolph Research DDM-2910 Densitymeter

- Bubble detection system using scanning video imaging and 10x magnification,
- Programmable method options,
- Measurements using less than 1 ml sample,
- LCD Touch-screen,
- Measurement Range: 0-3 g/cm³,
- Temperature Range: 0°C - 95°C,
- Resolution: 0.0001 g/cm³,
- Windows 7 operating system and 21 CFR Part 11 (User authorization) option.

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Application Areas of Density Meters

- Food Industry (Density measurements in alcoholic beverages),
- Academia (Research Laboratories),
- Chemical Industry (Analysis of chemicals),
- Pharmaceutical Industry (Raw material analysis),
- Petroleum Industry (Analysis of petrochemical samples),
- Agriculture Industry (Measurements of fertilizers -UAN- and agricultural chemicals),
- Cosmetics Industry (Measurement of fragrances).

Concentration and Solid Matter in Cosmetics

When light travels through a material of one density to one of different density, the direction change is called refraction, and the ratio of the speed of light velocity in a vacuum to the speed of light through the substance is called the refractive index. Glycerine, ethyl alcohol, acetone and other pure fluids have their own specific refractive index values.

Refractive index is a physical property of a substance and it is used in raw materials analysis as well as for the refractive index and concentration measurements of liquid and gel end products and those products showing fluid properties for quality control purposes in the cosmetic industry.

Rudolph Research J257 Refractometer

- Measurement Range: 1.26 - 1.72 RI, 0-100 BRIX,
- Measurement Accuracy: ± 0.0001 RI, ± 0.1 BRIX,
- Temperature Control: 10°C - 100°C, 8 GB memory,
- Touch-screen,
- Windows 7 operating system and 21 CFR Part 11 (User authorization) option,
- Smart Measure (Feedback on sample quality, correct sample loading and cleaning),
- Trend Analysis (allows the user to save a measurement into a method history to track long term stability of the user's unique methods).

Application Areas of Refractometers

- Cosmetic Industry (Measurement of concentration and solid matter),
- Academia (Research Laboratories),
- Food Industry (Sugar industry, olive oil, fruit juice and honey manufacturers. Purity measurements, solid materials dissolved in water),
- Chemical Industry (Refractive index measurements),
- Pharmaceuticals Industry (Quality control of raw materials and end products),
- Petroleum Industry (Oil ratios and concentration),
- Agriculture Industry (Solid matter measurements in fluids).