Ozan Halisçelik

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Assistant Service Manager

Ant Teknik Cihazlar

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Various applications such as flavours, preservatives, hyaluronic acid and SPF (sun protection factor) in cosmetics and perfumery products are of concern due to quality control, product development and R&D studies and regulatory purposes. Several analytical techniques such as HPLC, GC, GCMS, LCMSMS, UV-VIS-NIR are used to meet the demands from cosmetics industry.

Hyaluronic acid is an acidic mucopolysaccharide in which D-N-acetylglucosamine and D-gluconic acid are alternately bonded, and is known to be distributed throughout the body in tissues such as the skin, tendons, muscles, cartilage, blood vessels and brain; with a molecular weight in excess of 1,000,000. Since hyaluronic acid exists naturally in living organisms, its high biocompatibility is of interest, and it has been used in applications such as surgical procedures. Since hyaluronic acid gradually degrades in vivo, there is a growing need for hyaluronic acids with higher molecular weights, which remain effective for longer periods. Hyaluronic acid also has a high moisturizing effect, and thus is used as a cosmetic additive. Pharmaceutical formulation impurities and production step fractions of hyaluronic acid are controlled by LCMSMS ESI. Here, we would like to introduce examples of the use of Size Exclusion Chromatography (SEC) for the analysis of hyaluronic acid.

Analysis of Standard Samples

Fig.1 shows the structural formula of hyaluronic acid. When using size-exclusion chromatography for the analysis of macromolecules that contain ionic dissociation groups, these ionic groups within the molecule will repel each other, causing the molecule to expand. These ionic interactions can be suppressed by adding a salt in the mobile phase. In particular, in the case of a hydrophilic molecule such as polysaccharides, the existence of any sort of ionic groups will frequently cause pattern deformation and unstable measurement results. In this example, sodium sulfate was used as the salt. The analytical conditions are shown in Table 1.

Fig. 2 shows an example of the analysis of sodium hyaluronate derived from rooster comb. The sample is dissolved to give a 0.1% concentration in the 50mM aqueous solution of sodium sulfate used as the mobile phase.

Fig. 3 shows an example of the analysis of sodium hyaluronate derived from human umbilical cord. In the same manner as with Fig.2, the sample is dissolved to give a 0.1% concentration in 50mM aqueous solution of sodium sulfate.

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


Analysis of Cosmetics and Soft Drinks

We analyzed some cosmetics and soft drinks that claim to contain hyaluronic acid. The Analytical conditions are shown in Table 1. Chromatograms are shown for a soft drink in Fig.4, cosmetic A in Fig.5, and cosmetic B in Fig.6. Cosmetic A was analyzed after being diluted 10-fold in a 50mM aqueous sodium sulfate solution used as the mobile phase. The soft drink and cosmetic B were analyzed without dilution, but after filtration through a 0.45µm pore filter. The peaks between 11 and 17minutes are presumed to be hyaluronic acid.

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