Takayuki Yoshimori, Ph.D

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
In bio-pharmaceutical products, immunogenicity associated with an aggregation phenomenon of the protein has been attracting a lot of attention as a serious problem during clinical trials and after launch.

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Working range of Analytical Methods used for evaluating Protein aggregate/Insoluble particulate matters including Subvisible particles (SVP) are shown in Figure 1. As shown in this figure, Size Exclusion Chromatography (SEC) has been used to evaluate protein oligomers/ aggregate less than size of 0.1-0.2 μm . SEC method has many advantages over other methods such as its simple experimental protocol and high reproducibility. However, due to the presence of many problems (Decline of recovery ratio caused by trapping on the column top and absorption on the matrix & surface of the column, dissociation of protein oligomers/ aggregates by dilution of sample by mobile solution during analysis and its high pressure, etc.) in the quantitative evaluation of protein oligomers/ aggregates has received pointed out with respect to its legitimacy by research developers and FDA¹.

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
Figure 1. Working range of Analytical Methods used for evaluating Protein aggregate/Insoluble particulate matters including Subvisible Particles (SVP)

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For the evaluation of insoluble particulate matters over 10 um and 25 um, Light Extinction/ Obscuration analysis (Analytical device: HIAC) has been used in worldwide.² [Figure 2].

Dave Thomas et al, Particle Characterization Using Flow Microscopy, 2008 BioProcess International Conference & Exhibition


Şekil 2. Light Extinction/Obscuration Metodu Çalışma Prensipli ve Özellikleri
Figure 2. Principle and Characteristics of Light Extinction/Obscuration Method

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Although Light Extinction/ Obscuration method has been used for evaluating insoluble particulate matters over 10 um and 25 um in biopharmaceuticals for long time, there is serious problems for evaluating protein aggregates in biopharmaceuticals.

- As shown in Figure 3, since the detection sensitivity is proportional to the area and density of the shadow caused when light is applied to the particles, protein aggregates [Figure 3a] and silicone oil droplet [Figure 3b] are difficult to detect because those particles transmit light.
- Air bubble is highly sensitive for Light Extinction/ Obscuration method because air bubbles block the light as shown in Figure 3c.
- It is impossible to know what type of particles have been detected.

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

Figure 3. Transmission microscopic images of protein insoluble particle, Silicon oil droplet, and air bubble using microfluidic cell

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Indeed, in the measurement by Light Extinction/ Obscuration method, protein aggregates have been evaluated as smaller size compared to the actual particle size [Although depending on the thickness of protein aggregates, usually less than one-tenth of their actual size: data not shown]. To evaluate insoluble particulate matters more accurately, Microflow digital imaging method which has similar and even wider working range has been used as supporting/ alternative method of Light Extinction/ Obscuration method. Principle of Microflow digital imaging method is simple. Just take picture of all particles and evaluate the size of particle from the actual picture [Figure 4a & 4b]. Therefore, unlike Light Extinction/ Obscuration method, Microflow digital imaging method is able to evaluate mostly transparent protein aggregates so as silicon oil droplet. Replacement of Extinction/ Obscuration method to Microflow digital imaging method for evaluating insoluble particulate matters are very important to protect patients from large protein aggregates exist in biopharmaceuticals.

Dave Thomas et al, Particle Characterization Using Flow Microscopy, 2008 BioProcess International Conference & Exhibition

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Figure 4. Principle and Characteristics of Microflow Digital Imaging (4a) and images taken by Microflow digital imaging method (4b)

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
As indicated in Figure 1, protein oligomers/ aggregate less than 0.1-0.2 μm, insoluble particulate matters over 10 μm & 25 μm, and particles able to be detected by visual inspection are evaluated/ regulated. It means particles between 0.1-0.2 μm and 10 μm has not been evaluated/ regulated at present. Therefore, this size of particles is classified/ defined as Subvisible particles (SVP)³. Recently, the particles of this size exhibits a high immunogenicity has been confirmed [Figure 5]. One of the reason why SVP has not been evaluated is lack of analytical methods. Archimedes using Resonant Mass measurement, NanoSight using nanoparticle tracking analysis, and Aggregates Sizer

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using Laser Diffraction/ Scattering method are the devices specially made for evaluating SVP [Figure 1].

How does the Immune System Respond To Particulates $\leq 10\mu\text{m}$

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By the way, it was suggested that the definition of SVP may change from the current particle size range of 0.2 – 10 µm to particle size range of 1 – 100 µm at "2nd Progress and Challenges in Protein Particles and Immunogenicity of Biotherapeutics 2015" [Figure 1]. Since SVP analyzing devices such as Archimedes and NanoSight are not capable of evaluating new SVP range, and use of orthogonal methods are recommended by FDA⁴, use of Microflow Digital imaging method and Aggregates Sizer, which has quantitatively comparable measuring ability of same particles [Figure 6], are highly recommendable for evaluation of SVP with new SVP range.

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
Figure 6. Comparison of aggregates analysis results between Aggregate Sizer and MFI

Sample: IVIG solution (0.87 mg/mL, PBS pH 7.4)
Stress: Heat stressed for 15 min. at 70°C

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Kaynaklar:

1. Shire SJ. Assessment of Self-Association of Protein Therapeutics at High Concentration and its Impact on Viscosity. International Light Scattering Colloquium 2006.
2. Thomas D. Particle Characterization Using Flow Microscopy. 2008 BioProcess International Conference & Exhibition: Formulation Strategies for Protein Therapeutics
3. Carpenter JF. et al. 2008. Overlooking Subvisible Particles in Therapeutic Protein Products: Gaps That May Compromise Product Quality. J Pharm Sci, ODI: 10.1002/jps.21530.
4. Susan L. Kirshner. Regulatory Expectations for Analysis of Aggregates and Particles. 2014 Workshop on Protein Aggregation and Immunogenicity, 2014.July 17th.

Finally, quantitative evaluation methods of protein oligomers / aggregates, insoluble particles, and SVP in biopharmaceuticals have been described, and it is true that there are many problems in accuracy for evaluation methods using at present. As mentioned at the outset, immunogenicity occurs due to the aggregation phenomena of the biopharmaceuticals has become risk to the patient's safety. It can be said that Laser Diffraction / Scattering (Aggregates Sizer) technique for evaluating aggregates/SVP/insoluble particles in biopharmaceuticals is recommendable due to its wide measurement range (7 nm – 800 um).

Resources:

1. Shire SJ. Assessment of Self-Association of Protein Therapeutics at High Concentration and its Impact on Viscosity. International Light Scattering Colloquium 2006.
2. Thomas D. Particle Characterization Using Flow Microscopy. 2008 BioProcess International Conference & Exhibition: Formulation Strategies for Protein Therapeutics
3. Carpenter JF. et al. 2008. Overlooking Subvisible Particles in Therapeutic Protein Products: Gaps That May Compromise Product Quality. J Pharm Sci, ODI: 10.1002/jps.21530.
4. Susan L. Kirshner. Regulatory Expectations for Analysis of Aggregates and Particles. 2014 Workshop on Protein Aggregation and Immunogenicity, 2014.July 17th.