Sezen Çepel

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Analysis of Aroma Compounds and Fragrances by Using Flavour & Fragrance Natural & Synthetic Compounds (Ffnsc) Gcms Library

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Aroma compounds in food and beverages, and fragrances in cosmetics and perfumes are complex structures derived from natural and synthetic compounds.

GCMS is commonly used for flavor and fragrance analyses, and searching a library of GCMS spectra is convenient for identifying unknown components. However, typical GCMS libraries contain an insufficient number of registered flavor and fragrance components.

Additionally, the resemblance in mass spectra of some compounds causes difficulty in identification. The FFNSC library, created using Shimadzu's GCMS-QP2010 series by Prof. Mondello's group (Messina University, Italy), utilizes the GCMSsolution software function that allows users to sort candidate compounds by retention index, following a library search of mass spectra, making it ideally suited for flavor and fragrance analysis.

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Figure 1. GCMS-QP2010
Ultra Gas Chromatography
Mass Spectrometer

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
The measured spectrum (retention index 1090) and the spectra of compounds with a high score running a library search

Similarity Search by Using Retention Index

The data registered in FFNSC 2 contain the linear retention index (LRI) for every compound. The retention time, in addition to the mass spectrum, is important information for identification. GCMSsolution searches using both mass spectrum and retention index, significantly improving the reliability of identification. All data registered in FFNSC 2 libraries have been obtained using the Shimadzu GCMS-QP2010 series.

Shimadzu FFNSC (Flavour and Fragrance Natural and Synthetic Compounds Library) GC/MS Library FFNSC 2 has been expanded to include 3,000 mass spectra and retention index providing great support for flavor and fragrance analysis.

Filtering with Retention Index

Since flavor and fragrance components sometimes contain compounds having similar chemical structures, multiple compounds with similar mass spectra are listed as candidates when performing a library search using only the mass spectrum. Filtering with the retention index sorts the candidates by retention index, thereby providing highly accurate identification results.

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Searching both the MS spectrum and retention index enhances reliability of identification.

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Flavor/Fragrance Analysis by Fast-GC/MS

Fast-GC/MS is an effective way to improve laboratory productivity by shortening analysis cycle times. Fast-GC/MS methods require a system with a high-performance AFC controller that is compatible with narrow bore capillary columns and high-speed data acquisition technology. Shimadzu GCMS systems offer the high-performance levels necessary to satisfy these requirements. This application shows results from analyzing lavender oil using conventional methods and Fast-GC/MS methods.

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
Results

Fast-GC/MS provided separation patterns similar to results from conventional methods, but analysis times were 1/7 of conventional analysis times.

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GC/MS Quantitative Analysis of Allergens in Cosmetics Using OPTIC-4

Because an estimated 1 to 2 percent of the population has allergies, in 2003 the EU began regulating the use of perfumes and flavour materials that can cause allergic reactions. These substances are restricted to 0.01 % for

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cosmetic products that are rinsed away (such as shampoo) and 0.001 % for products that are left applied (such as hand creams). Due to these regulations, manufacturers now must quantitatively analyze the concentration of such allergens.

Experiment Sample

Standard samples of allergen mixtures and cosmetic cream diluted in solution by ten times were analized.

Analysis Conditions

Instruments

Inlet : OPTIC-4 (ATAS GL International BV)
Liner : L100011, DMI liner (ATAS GL International BV)
GC-MS : GCMS-QP2010 Ultra (Shimadzu)
Autosampler : AOC-5000 (Shimadzu)
LINEX (ATAS GL International BV)
Column : InertCap 5 MS/Sil (30 m x 0.32 mm I.D., df=0.5 um) (GL Science)

[Inlet]

Injection mode : Split
Balancing time : 5 sec,
Analysis time : 29.5 min
Vaporization chamber temperature : 35 C -> (5 C/min) -> 250 C (5min)
Carrier gas : Helium
Control mode : Flow control
Transfer time : 3 min,
Column flow rate : 1.0 ml/min,
Split flow : 100 ml/min.

[MS]

Interface temperature : 250°C
Ion source temperature : 200°C
Solvent elution time : 3 min,
Data sampling time : 3 – 39.5 min.
Measurement mode : Scan
Mass range : m/z 35-350
Detector voltage : Relative value mode


[GC]

Column oven temperature : 35°C (3 min) -> (10°C/min) -> 250°C (5 min).

Results

Standard Sample Measurement

Figures 3, 4 and 5 show the results of measuring 3 standard samples.

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
Figure 3. Measurement Result of Allergen Standard Sample 1

1. Benzyl Alcohol, 2. Linalool, 3. Citronellol, 4. Geraniol, 5. Anisyl Alcohol, 6. Cinnamyl alcohol, 7. Eugenol, 8. Isoeugenol, 9. Amylcinnamyl alcohol, 10. Farnesol, 11. Farnesol, 12. Cinnamaldehyde, .alpha. -hexyl -

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
Figure 4. Measurement Result of Allergen Standard Sample 2

1. Citral, 2. Citral, 3. Citronellal, hydroxy -, 4. Alpha.Isomethyl ionone, 5. Linal, 6. Cinnamaldehyde, .alpha. -pentyl - 7. Amylcinnamyl alcohol, 8. Lyril (21.111), 9. Lyril (21.238), 10. Cinnamaldehyde, .alpha. - hexyl

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Figure 5. Measurement Result of Allergen Standard Sample 3


1. Limonene, 2. Metil 2 - octynoate, 3. Coumarin, 4. Benzil Benzoate, 5. Benzil salicylate, 6. Benzil cinnamate

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Actual Sample Measurement

Cosmetic cream measurement results are shown in Figure 6.

Bileşik Adı Compound Name	Alıkonma Zamanı (dk) Retention Time (min)	m/z	Alan Area	RSD'ler RSDs	Miktar (µg/ml) Amount (µg/ml)
D-Limonene	12.241	68.00	2525392	%8	30.80
Citral	15.639	41.00	725549	%7	12.92
Geraniol	15.683	69.00	696831	%9	7.02
Citral	15.954	69.00	282420	%5	6.13
Cinnamyl Alcohol	16.634	92.00	51418	%5	1.07
Eugenol	17.246	164.00	101188.00	%6	1.60
.alpha. Isomethyl ionene	18.904	135.00	1563043	%9	13.05
Lilial	19.564	189.00	4009084	%5	40.10
Cinnamaldehyhde, .alpha.-pentyl-	21.000	129.00	370996	%6	7.09
Cinnamaldehyhde, .alpha.-hexyl-	22.145	129.00	3144191	%1	63.36
Benzyl salicylate	23.648	91.00	13512534	%3	67.06

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Figure 6. Cosmetic Cream Measurement Results

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Referanslar

- Shimadzu Uygulama Notu, GCMS No. 32
- Shimadzu Uygulama Notu, GCMS No. 37

The total ion chromatogram and area and repeatability accuracy of detected components are indicated above.

Measuring the allergens in cosmetic cream resulted in repeatability and % RSD values of 10 % or less, which indicated that the given method can be used to quantitatively analyze allergens.

References

- Shimadzu Application Note, GCMS No. 32
- Shimadzu Application Note, GCMS No. 37