Makale / Article

Zafer Bayram

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The Importance of Gas Chromatography In The Impurity Testing Of Process Gases Used In The Pharmaceutical Industry

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In the pharmaceutical industry, gas chromatography technique has its routine fields of application like analysis of volatile and semi-volatile compounds and residue solvent quantitation. Recently, this technique has gained even more popularity for another application: impurity analysis of high purity pharmaceutical process gases.

Since finished pharmaceutical products have direct effects on public health and environment, the pharmaceutical industry manufacturing process needs the utmost care and precision. Even very minor faults in any production stage may lead to severe adverse outcomes. Therefore, all stages of pharmaceutical processing including raw material income, packaging and transportation have to be controlled with the utmost care.

In the pharmaceutical industry, it is necessary to use several analytical techniques in order to check whether the raw materials and the products have the specified quality and product specifications. Gas chromatography (GC) is one of the most commonly used analytical techniques in the pharmaceutical laboratories. GC and GCMS are used majorly for the analysis of volatile and semi-volatile organic synthesis intermediates; active substances in the final product formulation or for the analysis of organic volatile impurities (OVIs) present in the raw materials, finished products and packaging materials. Recently, gas chromatography has gained even more popularity for another routine application; impurity testing of process gases used in the pharmaceutical industry.

During the pharmaceutical manufacturing process, pure gases or gas mixtures such as nitrogen, hydrogen, ammonia, carbon dioxide, boron hydride etc. are being used for several purposes including chemical synthesis, lyophilization, extraction, sterilization, cryogenic milling, hydrogenation, providing an inert atmosphere for avoiding oxidation

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and keeping thermal stability etc. Nitrogen is the most commonly used process gas used during production, storage and transportation stages in the pharmaceutical industry. Being an inert gas, nitrogen is used for providing a dry and nonflammable environment for storage and production. Other fields of use of nitrogen gas is freeze-drying, aseptic packaging, transportation of the products at a specific temperature or as a gas propellant for aerosols in the form of nitrous oxide.

Process gases used in the pharmaceutical industry are not classified under medical products and therefore GMP compliance for the manufacturing of these gases is not obligatory. On the other hand, starting from the early 90s, Food and Drug Administration (FDA) has focused on the nitrogen used for inerting the final product and this raised the demands from the gas suppliers. Other process gases used during the earlier stages of the pharmaceutical manufacturing on the hand receive lower interest. Pharmacopeias are strictly monitored and followed by the authorities, active pharmaceutical ingredient (API) and pharmaceutical manufacturers. Therefore, the quality of the process gases are evaluated according to the pharmacopeia monographs. US Pharmacopeia (USP) and European Pharmacopeia (EP) are the pharmacopeias which are most commonly referred to.

Nitrogen		
Monograph	EP	USP
Name of the Substance	Nitrogen	Nitrogen
Reference	1247	7727-37-9
Chemical Formula	N ₂	N ₂
Definiton	Min. % 99.5 V/V N ₂	Min. % 99.0 V/V N ₂
Identification	Retention time of peak with Gas Chromatography	Gas Chromatography currently under review
Production		
Assay	≥ % 99.5 V/V N ₂	≤ 1.0% O ₂ , ≥ % 99.0 V/V N ₂
Analytical Method	Gas Chromatography	Gas Chromatography
Impurities		
CO	≤ 5 ppm V/V	≤ 0.001 % V/V
CO ₂	≤ 300 ppm V/V	-
O ₂	≤ 50 ppm V/V	≤ 1.0 %
H ₂ O	≤ 67 ppm V/V	-

Table 1. Pharmacopeia monographs defined for nitrogen gas


Recently, impurity analysis for process gases are being realized by the pharmaceutical manufacturer besides the gas supplier.

GC-TCD (Thermal Ionization Detector) is used for this purpose. TCD can successfully be used in the analysis of substances like nitrogen, CO₂, formaldehyde, formic acid etc. which cannot be analysed by using FID. Helium or hydrogen are generally used as the carrier gas.

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Molecular sieve PLOT capillary columns are usually used as the column for this application. These columns are specially designed for efficient separation of gases like Ar/O₂, N₂, CH₄, C₂H₆ and CO. Special coating and deactivation procedures ensure chromatographic efficiency.

Figure 1 demonstrates an example of gas analysis using GC-TCD and Molecular Sieve 5A PLOT column.



Figure 1. Gas analysis application using GC-TCD and Moleküler Sieve 5A Plot Column

Şekil 3. Restek Molecular Sieve PLOT Kapiler Gaz Kromatografi Kolonu
Figure 3. Restek Molecular Sieve PLOT Capillary GC Column