Murat Onul

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
### Analysis of Benzalkonium Chloride in Textiles Using Lc/Ms/Ms

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*Benzalkonium chloride, also known as BZK, alkyl-dimethylbenzylammonium chloride and ADBAC, is a cationic surface-acting agent belonging to the quaternary ammonium group. It has three main areas of use; as a biocide, a cationic surfactant, and phase transfer agent in the chemical industry. The chemical is a nitrogenous mixture of alkylbenzyl-dimethylammonium chlorides of various even-numbered[1] alkyl chain lengths.*

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*Benzalkonium chloride has trade names as BC50, BC80, BAC50 and BAC80.*

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#### Some Areas of Use

- *In textile industry; as protection agent for natural fibers against mold, fungi and moth.*
- *In textile and leather (shoes) industries; added during the last steps of the process as an antibacterial agent.*
- *In leather, textile and paper industries; during wetting and as retarder during dyeing process of acrylic fibres using cationic dyes.*
- *In leather industry; for retardation of mold formation and for easier softening, wetting and dyeing of the leather.*
- *In pharmaceutical industry; as skin antiseptics and sanitizers.*
- *In cosmetics; as new generation hand cleansing disinfectants and anti-microbial agents.*
- *In swimming pools; as disinfectants (water and effluent treatment) and as algacides.*

*%10 (or more) solution of Benzalkonium chloride is toxic*

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for human and may cause irritation in skin and mucosa. Scientific studies have proven the genotoxic effect [2] of benzalkonium chloride and have shown that it might cause allergic reactions. It is important that the content of benzalkonium chloride is verified.

In this article, Shimadzu LCMS8030 Triple Quadrupole system has been used for the fast and simple quantitative analysis of benzalkonium chloride in the textile industry.

Şekil 1. LCMS-8030 Triple Quadrupole Kütle Spektrometre Sistemi  
Figure 1. LCMS-8030 Triple Quadrupole Mass Spectrometer

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**Table 4. MRM values and collision energies (CE) for compounds.**

| Compound Name | MRM Transitions  | CE  |
|---------------|------------------|-----|
| BAC8          | 248.2/91.1/156.4 | -15 |
| BAC10         | 276.2/91.1/184.4 | -12 |
| BAC12         | 304.2/91.1/184.4 | -15 |
| BAC14         | 332.2/91.1/240.3 | -10 |
| BAC16         | 360.2/91.1/268.4 | -12 |
| BAC18         | 388.2/91.1/296.4 | -15 |

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| No | Bileşik İsimleri | Q1    | Q3    | RT (Alıkonma Zamanları) | %RSD (20pp, n=6) |
|----|------------------|-------|-------|-------------------------|------------------|
| 1  | BAC 8            | 248.1 | 90.9  | 1.98                    | 1.12             |
| 2  | BAC 12           | 304.2 | 91.1  | 4.13                    | 2.31             |
| 3  | BAC 14           | 332.2 | 91.1  | 6.81                    | 2.46             |
| 4  | BAC10            | 276.1 | 91.0  | 4.51                    | 0.75             |
| 5  | BAC 16           | 360.3 | 91.1  | 8.10                    | 1.11             |
| 6  | BAC 18           | 388.3 | 129.1 | 24.7                    | 0.81             |

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

**Table 5. MRM values, retention times and relative standard deviation percentages for the compounds at 0.1 ppb concentration.**

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#### Referanslar / References:

1. U.S. Environmental Protection Agency: Reregistration Eligibility Decision for Alkyl Dimethyl Benzyl Ammonium Chloride (ADBAC).
2. Benzalkonium chloride (BAC) and dimethyldioctadecyl-ammonium bromide (DDAB), two common quaternary ammonium compounds, cause genotoxic effects in mammalian and plant cells at environmentally relevant concentrations, F. Ferk, M. Mišík<sup>1</sup>, C. Hoelzl, M. Uhl<sup>2</sup>, M. Fuerhacker<sup>3</sup>, B. Grillitsch<sup>4</sup>, W. Parzefall, A. Nersesyan, K. Mičeta<sup>1</sup>, T. Grummt<sup>5</sup>, V. Ehrlich and S. Knasmüller.\*