

Energy Dispersive X-ray Fluorescence Spectrometer

# ALTRACE



# ALTRACE™

ENERGY DISPERSIVE X-RAY FLUORESCENCE SPECTROMETER

## Pushing the Boundaries of Detection

Detect trace elements with ease.

A combination of optical system design and Shimadzu's proprietary high-speed signal processing technology allows ALTRACE to reach new heights in terms of sensitivity.

## Unparalleled Sensitivity

- Implement batch elemental analysis from the ppm to percent levels.
- Screen from sub ppm to percent levels.
- Convenient analysis of liquid and powder samples.

## Escape from Complicated Pretreatment

- No need for complicated chemical pretreatment before analysis.
- Direct and simple screening, and also precise.
- Lower cost and easier operation compared to methods that incorporate pretreatment.

## High Throughput and Efficiency

- Automatic consecutive analyses of up to 48 samples
- User-friendly drawer style tray
- Support for interrupted analysis during measurements

 SHIMADZU

ALTRACE ENERGY DISPERSIVE X-RAY FLUORESCENCE SPECTROMETER

SHUTTER OPEN

X-RAYS ON

ERROR

# Unparalleled Sensitivity


With Shimadzu's combination of optical system design and proprietary high-speed signal processing, ALTRACE reaches the highest detection sensitivity compared to Shimadzu's general-purpose systems (EDX-7200).

## Unparalleled Sensitivity that Goes Beyond Typical Measurement Needs


The combination of optical system that optimally arranges high-power X-ray tube and high-sensitivity detector, and high-speed signal processing has achieved high detection efficiency.

In the analysis of aqueous solution samples, the detection limit for all elements have been improved compared to the conventional equipment (EDX-7200).


It can be used effectively for the analysis of low concentration samples and the automatic analysis of multiple samples requiring shorter measurement time.




Comparison of detection limit between conventional equipment (EDX-7200) and ALTRACE




EDX-7200  
(600 seconds measurement time)  
Pb 5ppm



ALTRACE  
(600 seconds measurement time)  
Pb 1ppm




EDX-7200  
(600 seconds measurement time)  
Cd 5ppm



ALTRACE  
(600 seconds measurement time)  
Cd 1ppm





Implements batch element analysis across a wide range from the sub ppm level to the % level



Automatic Filter Exchange Improves Sensitivity

Background removal is facilitated by a primary X-ray filter, thereby improving S/N ratio. For trace-level analysis, where background is not negligible, filters are particularly advantageous. ALTRACE is equipped with 6 primary filters (8 types, including open and attenuator), which can be exchanged autonomously by the software.

	Effective Energy (keV)	Example of Applicable Elements
Filter #1	14 to 38	Mo,Rh,Pd,Ag,Cd,Sn,Sb
Filter #2	2 to 4	S,Cl
Filter #3	5 to 10	Cr,Mn,Fe,Co,Ni
Filter #4	9 to 10	Hg
Filter #5	9 to 14	As,Br,Zr,Hg,Pb,Bi
Filter #6	4 to 5	Ti,V




# Escape from Complicated Pretreatment

Liquids and powders can be directly measured by simply placing them as is in the sample cell. EDX is a non-destructive technique; therefore, the same sample can be measured for both simple screening and precision analysis. To date,

previous EDX systems could not measure concentrations on the order of 0.1 ppm. With ALTRACE, this level is now detectable. It is now possible to obtain sub-ppm level detection, without the need for chemical pretreatment.

## Simple Screening in Combination with Precision Analysis

When performing quantitative analysis by specifying elements, once a calibration curve is created, there is no need to recreate the calibration curve for each measurement. In addition, ALTRACE is suitable for simple screening analysis, since qualitative and quantitative analysis using the fundamental-parameter (FP) method can be performed without specifying elements.



## Solids, Powders and Liquids Can Be Measured As Is

With ALTRACE, samples are placed into a dedicated sample cell. The sample is supported by a transparent X-ray film that is suitable for the sample.

### Powdered Samples and Pellet Samples

Powder and pellet samples are placed in the sample cell and supported with polypropylene film.



Powder

Pellet

### Liquid Samples (Water Soluble and Oleaginous)

Aqueous solutions are poured into the sample cell and supported on a polypropylene film. Other liquids, such as organic solvents, can be measured by using an appropriate film.



Aqueous solution


Organic solvent



## High Throughput and Efficiency

### Automatic Consecutive Analyses of Up to 48 Samples


710mm wide main body is equipped with a multi-sample changer. The tray drawer system allows for safe sample replacement without having to access the inside of the instrument.



**48**  
samples

### Analysis can be Paused and Reassigned, Increasing Flexibility

During analysis, ALTRACE can be paused, and new samples can be added and the analysis schedule adjusted accordingly. As such, there is no need to wait for analysis to be completed to start sample preparation.



# Ability to Measure a Diverse Array of Matrices Makes ALTRACE Suitable for a Variety of Applications



## Food Products Field

Cacao powder, rice, powdered milk, etc

Calibration Curve  
Method

Cadmium, arsenic, lead, and mercury are highly regulated in food products due to their inherent toxicity. A maximum allowable value of less than 1 ppm is specified. ALTRACE can detect these levels, and as low as 0.1 mg/kg, in food products.

### Analysis Results for Standard Sample NMIJ7502-a



n	Cd
1	0.508
2	0.599
3	0.597
4	0.474
5	0.557
6	0.549
7	0.469
8	0.502
9	0.574
10	0.485
Mean	0.532
Standard deviation	0.050
Standard value	0.548
Coefficient variation	9.40

NMIJ CRM 7502-a  
Results of 10 replicates [mg/kg]  
(300 seconds measurement time)




## Chemical Products and Soil

Raw materials for chemical products (liquids and powders)

FP Qualitative and  
Quantitative Analysis

Energy dispersive X-ray fluorescence spectrometers are used in fields related to the environment. Powdered and granular samples can be placed in a sample cell and analyzed as is, without dissolving them. Soil contains a variety of elements. After performing a qualitative analysis of such elements, a quantitative analysis is performed using the FP method.

### Analysis Results for Soil Samples



Target Elements	Analysis Results	Target Elements	Analysis Results
SiO <sub>2</sub>	59.793 %	ZnO	0.015 %
Al <sub>2</sub> O <sub>3</sub>	22.177 %	Rb <sub>2</sub> O	0.008 %
Fe <sub>2</sub> O <sub>3</sub>	8.074 %	CuO	0.005 %
K <sub>2</sub> O	1.472 %	SrO	0.005 %
TiO <sub>2</sub>	0.832 %	As <sub>2</sub> O <sub>3</sub>	0.003 %
MgO	0.715 %	NiO	0.003 %
CaO	0.417 %	Ga <sub>2</sub> O <sub>3</sub>	0.003 %
P <sub>2</sub> O <sub>5</sub>	0.268 %	Y <sub>2</sub> O <sub>3</sub>	0.003 %
SO <sub>3</sub>	0.079 %	PbO	0.002 %
MnO	0.063 %	NbO	0.001 %
V <sub>2</sub> O <sub>5</sub>	0.051 %	I	0.001 %
ZrO <sub>2</sub>	0.030 %	Br	0.001 %
Cr <sub>2</sub> O <sub>3</sub>	0.026 %		

Results of the quantitative analysis



# Consumables and Options

## Polyester Film

500pcs/set

202-86501-59

For supporting samples. Polyester is recommended for organic solutions, e.g. oils.

## Polypropylene Film

500pcs/set

212-14030-91

Film for supporting samples. Polypropylene is suggested for liquid powder, and pellet samples.

## Sample Cell Body

100 pcs/set, I.D. 25mm, volume 10mL

212-14043-01

Polypropylene



## Sample Cell Ring

100 pcs/set

212-14043-03

For fastening the sample film.



## Sample Cell Lid (for Air)

100pcs/set

212-14043-02

Lid of sample cell Dedicated part for auto-sampler



## Sample Cell Lid (for Vacuum)

100 pcs/set

212-14043-04

Lid of sample cell for use in vacuum



## Sample Plate (for Air)

212-13834-91

Additional sample plates (for air).  
Four sample plates (for air) are included as standard with the instrument.



## Sample Plate (for Vacuum)

212-13834-92

Additional sample plates (for vacuum).  
Four sample plates (for vacuum) are included with the vacuum measurement unit.



## Contamination-Proof Cup

10 pcs/set

212-14114-41

Protects the detector and X-ray tube from contamination.



## Vacuum Measurement Unit

212-14200-41 (for 100V)

212-14200-42 (for 220-240V)


For high-sensitivity measurement of light elements. Space behind the unit will have to be available for rotary pump and vacuum control switch box.

# Analysis Functions

## Calibration Curve Method

In the calibration curve method, standard samples are measured, and the relationship between standard concentration and X-ray intensity is manifested in a curve. The curve can be used to quantify samples of unknown concentration.

Standard samples should be composed of a similar matrix as the unknown. In addition, a curve must be made for each target element. Calibration curve methods can provide high accuracy and quantitative data. Correction methods, such as absorption/excitation correction, overlapping correction, and other methods, are supported.




## FP Method

The FP method determines the elemental composition using a theoretical intensity calculation. Useful for quantitative methods where it is difficult to procure standards. For samples where the main components are C, H, and O, a balance (residue) setting is required, which can be determined by the user beforehand or theoretically using the shape of the scattered X-ray profile.

## Background FP Method


The background FP method uses the intensity of scattered X-rays (background) to further correct for matrix effects.

[cps/uA] Ti-U



## Principle Behind the Production of X-Ray Fluorescence


When a sample is irradiated by X-rays, characteristic secondary X-rays are emitted by atoms within the sample. These secondary X-rays are referred to as fluorescent X-rays. Each element has its own unique set of characteristic X-rays with their own energy and wavelength. As such, an energy dispersive X-ray spectrometer can be used to do both qualitative and quantitative analysis by measuring the intensity of X-rays at a given energy.



## Main Specifications

<b>Measurement Specifications</b>	<b>Measurement principle</b>	X-ray fluorescence spectrometry
	<b>Measurement method</b>	Energy dispersive
	<b>Target samples</b>	Solids, liquids, powders
	<b>Measuring range</b>	Na to U
	<b>Sample cell size</b>	Dia. 25 × 20 mm
<b>X-Ray Generator</b>	<b>X-ray tube</b>	Rh target (100 W)
	<b>Voltage</b>	5 to 65 kV
	<b>Current</b>	4000 µA max.
	<b>Cooling method</b>	Air-cooled (with fan)
	<b>Analysis region</b>	Dia. 18mm
	<b>Primary filters</b>	Automatic exchange of 6 types (8 types when including open and the attenuator)
	<b>Leakage dose rate</b>	0.5 µSv/h max.
<b>Detector</b>	<b>Type</b>	Silicon drift detector (SDD)
	<b>Liquid nitrogen</b>	Not required (electronic cooling)
<b>Measurement Chamber</b>	<b>Measurement atmosphere</b>	Air, vacuum
	<b>Sample replacement</b>	48-sample autosampler
	<b>Protective function</b>	Contamination-proof cup
	<b>Sample rotation function</b>	Sample spinner included
<b>Data Processor</b>	<b>CPU</b>	Intel® Core™ i5-13500 or above
	<b>Memory</b>	16 GB min.
	<b>SSD</b>	500 GB min.
	<b>Optics</b>	Super Multi Driver
	<b>PC OS</b>	Windows® 11 Pro
<b>Software</b>	<b>Qualitative analysis</b>	Measurement/analysis software
	<b>Quantitative analysis</b>	Calibration curve method, matrix correction, FP method, thin-film FP method, background FP method
	<b>Utilities</b>	Automatic calibration function (energy calibration and intensity check)
	<b>Others</b>	Instrument status monitoring function, analysis results tabulation function
<b>Installation</b>	<b>Temperature</b>	12 to 30 °C (temperature fluctuation rate 2 °C/hour max., temperature fluctuation range: ±5 °C max.)
	<b>Relative humidity</b>	40 to 70 % (No condensation)
	<b>Power supply</b>	100 to 240 V ± 10 %, 50/60 Hz, 300 VA*
	<b>Dimensions</b>	W710 × D550 × H570 mm
	<b>Weight of main unit</b>	Approx. 110 kg

\* Options are not included.



The vacuum measurement unit (optional) consists of a rotary pump and a control switch box.  
Make sure that the rear and sides of the ALTRACE unit are at least 800 mm away from walls and other devices that are difficult to move.

ALTRACE is a trademark of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.  
Intel and intel Core are trademarks of Intel Corporation or its subsidiaries.  
Windows is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.